interTRAX™ Mobile ICS and EMS Software

User Guide Version 3 2011



Copyright © 2003 - 2011 by Salamander Technologies, Inc.1 Protected by US Patents 5,573,278; 5,596,652; 5,793,882; 6,761,312; 7,191,934

Contents

| C | ontents | 2 | , |
|---|--|------|---|
| | Disclaimer | 4 | ļ |
| 1 | Introduction to <i>interTRAX</i> TM Mobile Software | 5 | į |
| | Overview | | |
| | interTRAX TM Barcode & Smart Card Information | 6 |) |
| | interTRAX TM Barcode & Smart Card Encoding – Personnel | | |
| | interTRAX TM Barcode & Smart Card Encoding – Personnel Medical | 7 | 1 |
| | interTRAX TM Barcode & Smart Card Encoding – Evacuee (rapidTAG evac) or | | |
| | Patient (Mobile EMS) | | |
| | interTRAX TM Barcode & Smart Card Encoding – Equipment & Assignment | | |
| | interTRAX TM Barcode & Smart Card Encoding – Company | . 10 |) |
| 2 | interTRAX TM Mobile (ICS & EMS) | | |
| | Settings | | |
| | Initial Startup | | |
| | Servers Tab | | |
| | Add or Edit Server | | |
| | Set Server as Active | | |
| | Organization Tab | | |
| | Scanner Tab | | |
| | Time | | |
| | Basic Tab | | |
| | Advanced Tab | | |
| | Client tab | | |
| | Sync | | |
| | Incident Screen Functionality – ICS tab | | |
| | 1 Details | | |
| | Reassignment | | |
| | 2 Timer | | |
| | 3 PAR | | |
| | 4 Assignment | | |
| | 5 Incident Details | | |
| | 6 Menu | | |
| | Menu prior to an incident | | |
| | New Incident | | |
| | New Medical Incident | | |
| | Settings | | |
| | Time | | |
| | <i>Sync</i> | | |
| | About | | |
| | Menu during an Open Incident | | |
| | Details | | |
| | Find | | |
| | Add Companies | . 25 |) |

| Roster | 25 |
|--|----|
| 7 Soft Keyboard | 25 |
| 8 Sync | 25 |
| 9 OK | 25 |
| Roster | 26 |
| Add a New Company to the Roster - Manual | 26 |
| Enter Company Information | |
| Add Resources to a Company in the Roster – Manual | 28 |
| Add Resources to a Company in the Roster – Manual | 28 |
| Disband a Company | |
| Remove Resources from a Company in the Roster | 30 |
| Access the Roster during an Open Incident | |
| Using interTRAX TM MOBILE ICS | 33 |
| Start a New Incident | 33 |
| Add Resources to an Open Incident | 34 |
| Add a New Company to an Open Incident - Manual | 35 |
| Add a New Company to an Open Incident - Electronic | |
| Remove Resources from a Company in an Open Incident - Manual | 38 |
| Remove Resources from a Company in an Open Incident - Electronic | 39 |
| MOBILE ICS Functionality | 40 |
| Examples | |
| Open, Close and End an Incident | 42 |
| Open | |
| Close or end on incident | |
| Using interTRAX TM MOBILE (EMS) | |
| New Incident` | 43 |
| Detailed | 44 |
| Scan | 44 |
| Manual Entry | 44 |
| Rapid Entry | 45 |
| Patient Details | 46 |

Disclaimer

Salamander Technologies, Inc. has made its best efforts to make the *Salamander* system reliable and defect-free; however, any computerized system is subject to malfunction including operator error, undetected system "bugs", and power or communication outages. Salamander Technologies, Inc. cannot and does not warrant that the system will perform flawlessly at all times and under all conditions.

1 Introduction to *interTRAX*[™] **Mobile Software**

Overview

Welcome to the *interTRAX*TM Mobile User Guide. *interTRAX*TM Mobile software can be used to collect data from Salamander tags. The information that is collected can be sent to Command software while on site.

Salamander is setting the standard for resource accountability in the greater U.S. Fire Service and was designed with three objectives in mind:

- **Simplicity** Salamander uses nationally recognized incident management procedures. Use of the Mobile Unit [MU] is self-explanatory and intuitive. Mobile Units, often referred to as Pocket Personnel Computers, set the standard for simplicity in the industry. Most important, the PDF417 ID tags carry the person's database just aim and shoot the MU to capture and display this data (name, rank, personal ID number, department, qualifications, etc.) on the spot. Operators can be trained to use the system in a matter of minutes.
- **Survivability** There are three tiers of survivability in Salamander System so that even a complete breakdown in electronic equipment does not inhibit accountability:
 - If communications to Central Dispatch or the internet are lost, incidents can be run from the Command Module.
 - If the Command Module is inoperable, basic accountability can be maintained via the MU's.
 - In the unlikely event that the MU's are down too, ID cards can be used as a traditional dog tag system to allow manual accountability.
 - PDF417 barcodes represent the ultimate in a distributed database they survive even if the Host system is destroyed or inoperable.
- Scalability The *interTRAX*TM Suite was designed to scale across the entire fire service so that it would be the right size for the needs and budgets of any department. Small departments can start with an *interTRAX*TM Suite compatible dog tag system. Medium departments can purchase a basic *interTRAX*TM Suite system for only a few thousand dollars. Large departments can deploy multi-MU systems as part of a Wireless LAN controlled by a Command Module. You can barcode only your department or you can barcode all mutual aid departments in the surrounding area. If the incident expands beyond mutual aid, you can use a label printer to barcode new responders at staging. Tell your mutual aid partners about the system: The *interTRAX*TM Suite unites fire departments, mutual aid, and state & national responders with a common machine-readable format.

interTRAX™ Barcode & Smart Card Information

The *interTRAX*TM barcode is a two dimensional, high capacity barcode that follows the PDF417 data format. The *interTRAX*TM smart card is an ISO standard 14443 or ISO standard 15693 contactless RFID card. The use of the information contained within either medium is protected by United States of America and International patents.

Within the *interTRAX*TM suite the contents of the Smart Card and the barcode are identical.

The first time a Smart Card is encoded it is setup to hold either company or non-company information. After the initial setup of a company Smart Card it may be encoded with company information and a non-company Smart Card may be encoded with assignments, personnel, equipment and victims.

interTRAX™ Barcode & Smart Card Encoding – Personnel

Personnel Container

The personnel container holds the following information. The items in bold print may be displayed in human-readable text. The human-readable text that appears is dependent on the type of ID tag, label or hardcard. The different types are listed below with the human-readable text that may appear if it is held in the container. The information available in the container and in human-readable text is dependent on the information that has been entered into the person's record.

- 1. Organization Country
- 2. Organization State
- 3. Organization Type
- 4. Organization ID
- 5. Organization Name
- 6. Personnel ID
- 7. Last Name
- 8. First Name
- 9. Rank
- 10. Date of Birth
- 11. Qualifications

Human-Readable Text

Velcro Tag Label

- First Name
- Last Name
- Personnel ID

Personnel Label

- First Name
- Last Name
- Organization ID
- Personnel ID

• Organization Name

CR-80 PVC Card

- Organization Name
- First Name
- Last Name
- Rank
- Organization ID
- Personnel ID
- Qualifications

Laminated Manual ID Tag

- Organization Name
- First Name
- Last Name
- Rank
- Organization ID
- Personnel ID
- Qualifications

Avery Label

- First Name
- Last Name
- Rank
- Organization ID
- Personnel ID
- Organization Name

interTRAX™ Barcode & Smart Card Encoding - Personnel Medical

Personnel Medical Container

The personnel medical container holds the following information. The items in bold print may be displayed in human-readable text. The information available in the container and in human-readable text is dependent on the information that has been entered into the person's record.

- 1. Organization Country
- 2. Organization State
- 3. Organization Type
- 4. Organization ID
- 5. Personnel ID
- 6. Physician Name
- 7. Physician Phone
- 8. Insurance Carrier
- 9. Insurance Policy number
- 10. Emergency Contact
- 11. Emergency Contact Phone
- 12. Organ Donor
- 13. Resting Pulse
- 14. Blood Pressure
- 15. Blood Type

- 16. Gender
- 17. Religion
- 18. Medications
- 19. Allergies

NOTE: The data read from the medical container is not retained by the MU's.

Human Readable Text

Laminated Manual ID Tag

The human-readable medical information is printed on the inside of the tag and is protected inside a heat-sealed laminate pouch. It includes all of the highlighted information above.

interTRAX[™] Barcode & Smart Card Encoding – Evacuee (rapidTAG evac) or Patient (Mobile EMS)

Evacuee/Patient Container

The evacuee/patient container holds the following information. The items in bold print may be displayed in human-readable text depending on the format of the ID Tag. The information available in the container and in human-readable text is dependent on the information that has been entered into the person's record.

- 1. First Name
- 2. Last Name
- 3. ID
- 4. Date of Birth
- 5. Gender
- 6. Class (e.g., adult, caregiver)
- 7. Triage Priority
- 8. Address
- 9. Special Needs
- 10. Parent
- 11. Triage Time

interTRAX™ Barcode & Smart Card Encoding – Equipment & Assignment

Equipment Container

The equipment container holds the following information. The items in bold print are displayed in human-readable text. The human-readable text that appears is dependent on the type of ID tag, label or hardcard. The different types are listed below with the human-readable text that may appear if it is held in the container. The information available in the container and in human-readable text is dependent on the information that has been entered in the record for the piece of equipment.

- 1. Organization Country
- 2. Organization State

- 3. Organization Type
- 4. Organization ID
- 5. Organization Name
- 6. Equipment ID
- 7. Manufacturer Make
- 8. Model
- 9. Description
- 10. Manufacture Date
- 11. Equipment Type
- 12. Date in Service
- 13. Features (qualifications)

Equipment Label

- Description
- Organization ID
- Equipment ID
- Organization Name

CR-80 PVC Card

- Organization Name
- Equipment ID
- Description
- Manufacture Date
- Manufacturer Make
- Model

Avery Label

- Description
- Organization ID
- Equipment ID
- Organization Name

Assignment Container

The assignment container holds the following information. The item in bold print is displayed in human-readable text. The information available in the container and in human-readable text is dependent on the information that has been entered in the record for the assignment.

- 1. Organization Country
- 2. Organization State
- 3. Organization Type
- 4. Organization ID
- 5. Assignment

Assignment Label

Assignment

Avery Label

• Assignment

interTRAX™ Barcode & Smart Card Encoding – Company

Company Container

The company container holds the following information. The items in bold are required.

- 12. Organization Country
- 13. Organization State
- 14. Organization Type
- 15. Organization ID
- **16. Organization Name**
- 17. Company Name
- 18. Company ID
- 19. Kind
- 20. Designation
- 21. NIMS Type
- 22. Qualifications
- 23. Leader Flag
- 24. Personnel Collection
- 25. Equipment Collection'
- 26. Crews (not in barcode)

Default Human-Readable text

- Company Name
- Company ID
- Organization
- NIMS Type
- Company Leader Name
- Qualifications

2 interTRAX™ Mobile (ICS & EMS)

The *interTRAX*TM Mobile software is used to:

- collect data from *Salamander* tags, smart cards and Triage Tags
- manually enter information
- send the information collected to *interTRAX*TM

NOTE: This User Guide should be used in conjunction with the Mobile Hardware Setup Guide where information on mobile unit screen setup and network setup can be found.

Tap on



in the lower left hand corner of the Main Screen



The interTRAX suite MOBILE screen will appear



Initializing...

Settings

Initial Startup

Following the installation of MOBILE the Settings screen will appear. It contains three tabs. You must complete the information on these tabs prior to using all of the functionality available in the MOBILE software.

The Settings screen does not appear each time you start MOBILE. If you want to make changes use the popup menu in the lower left hand corner of the main screen.


To change this information after you have started an incident close - not end - an incident and access the settings using the popup menu. Then open the incident.


$Menu \rightarrow Settings$




Servers Tab


Add or Edit Server









When the test is complete the screen will display a note that the server was connected







Organization Tab

Enter

Name. Type in the name of the organization

Code. The Code is a unique identifier and must be entered correctly. Contact STI Tech Support to

obtain the code.


Type. Choose from the dropdown **Country.** Choose from the dropdown **State**. Choose from the dropdown

Tap on Save



Scanner Tab

Default settings



Time


$Menu \rightarrow Time$



Basic Tab


Enter


- 1. IP Address for Command
- 2. An Internet Time Server (e.g., time-a.nist.gov) only if device can attach to the internet.



Advanced Tab

Uncheck the Phase Correction boxes






Sync

$Menu \rightarrow Sync$



The server settings are made when the Mobile Software is initially setup and may be changed at any time.

The settings should be tested prior to syncing to ensure that the syncing process is not delayed during an incident.





Incident Screen Functionality – ICS tab

The Incident screen contains 9 important elements

1 Details MOBILE Timer 00:00:00 PAR 4 2 Timer Assignment On Scene 3 PAR COMPANDED (1) 4 Assignment On Scene (4) - 😽 Water Tender 3 (4) **5 Incident Detail** 👸 DeRosa, Nick 6 Menu 👸 Hodge, Helen 7 Soft Keyboard 👸 Hogg, Jean 👸 Wentworth, Muriel 8 Sync MFD WT2 [WT2-03] 9 OK

1 Details




Highlight the company or the resource and tap on the **Details** icon to review the history for a company or the history and qualifications for a resource.



Reassignment

A resource may be reassigned from this screen by using the assignment dropdown or typing the assignment into the Assignment box.

Type in an assignment or choose the assignment from the dropdown



The resource will be moved to the new Assignment

2 Timer

The Timer is automatically set for 20:00 minutes

Minutes and Seconds

The minutes and seconds can be adjusted using the arrows or the soft keyboard

Auto Restart

Check the Repeat box if you want the time to restart automatically after it has timed out

Start and Stop

Tap on **Start** and **Stop** as appropriate to your needs

Time Expired


When the Timer expires the **PAR** screen will appear

Check the boxes appropriate for the information received when calling for PAR

When you check the box for a Company, the Responders are also checked

Checking all of the boxes for the Responders in a Company will automatically place a check in the Company box.






3 PAR

The PAR (Personnel Accountability Report) is the total number of personnel scanned on the scene. This number automatically increments up and down as personnel are scanned into Assignments and off of the scene.

4 Assignment

The Assignment dropdown includes all of the Assignments that have been scanned or manually entered into the MU. To manually enter an assignment, type it in the assignment box and click on the enter key or scan a resource and the assignment will be placed on scene with the resource attached. You may also add assignments on the Resource and Company Detail screen.



5 Incident Details

The Incident Details are presented in an incident command tree format

Assignments are not present in the tree when they are empty.



All Assignments that are added to an incident are placed subordinate to On Scene.

Each Assignment has a PAR count, shown in (), that represents all subordinates to that assignment.

On Scene () PAR represents all of the personnel on scene and matches the PAR on the upper right side.

Tap on in to expand the Assignment and view the Resources within

6 Menu

Menu prior to an incident

New Incident

See the New Incident Section

New Medical Incident

See the New Medical Incident Section

Settings

See the Options Section

Time

See the Time Section



Sync

The server settings are covered when the Mobile Software is initially setup and may be changed at any time. The settings should be tested prior to syncing to ensure that the syncing process is not delayed during an incident.

About



Menu during an Open Incident




Details



Details - Incident

Make changes in the incident information



Details - PAR

This is the same PAR check that is available when the timer is expired and is explained in the Timer section above.

Details - Medical

Read medical information from a smart card. NOTE: This screen must be open prior to reading medical information from a smart card.



Find



<u>Find – Responders, Equipment and Companies</u>

Check the Qualifications

Tap on **Find**

Only Resources or Companies with all of the Qualifications checked will be found

Add Companies

See the Add a New Company to an Open Incident section

Roster

See the Roster section

7 Soft Keyboard

Tap on the Soft Keyboard icon and at will appear at the bottom of the screen.

8 Sync

Tap on Sync and a series of screens will provide the status of the sync.






9 OK

Tap on **OK** in the bottom right hand corner and the End Incident box will appear.

Tap on **Yes** to end the incident.

Tap on **No** to close the incident, keep it IN PROGRESS and maintain PAR.




Tap on **Cancel** if you have reached this screen in error and do not want to end or leave Mobile.

Roster

Prior to starting an incident you may want to consider adding companies to the Roster. When you open an incident you will have access to the information in the Roster and may easily add it to the incident.

Add a New Company to the Roster - Manual

Prior to the start of an incident you may add companies to the Roster. This may be done manually using the menu or electronically by scanning barcodes.




Enter Company Information

Name. Type in a company name which is required. The maximum number of characters is 50.


ID. Type in a company ID which is required. The maximum number of characters is 50.

NIMS Type

Choose the appropriate information from the dropdown




The company will appear




Add Resources to a Company in the Roster - Manual


Highlight a Company


Scan in a resource




Add Resources to a Company in the Roster - Manual









Set Leader

Highlight the resource

Tap on the **Leader** box



The Leader will appear at the top of the list


Checked In

Check the boxes for the company and any resource that you wish to have the status of On Duty.





Tap on OK


The information will be stored in the Roster and will be available for use during an open incident.



Remove Resources from a Company in the Roster









Select Companies and/or Resources

From the Roster you may choose companies with or without the resources and any of the resources that are listed. The companies and resources that are checked will be added to the active incident assignment.






Tap on **OK**


Company with Resources



Select All or Select None



When you **Select All** the companies and all of the resources will be checked and added to the active incident assignment.




Using interTRAX™ MOBILE ICS

Salamander barcodes including those found on resourceMGR ID tags, Equipment ID tags rapidTAG badges and rapidTAG evacTM wristbands as well as triage tags may be scanned and used in MOBILE ICS.


Start a New Incident

The Main Screen will appear when the Settings have been entered during the initial setup and each time you start MOBILE.




The New Incident screen will appear

Name. Enter an incident name or leave it blank. **Number.** Enter an incident number or leave it blank **Type.** Choose from the dropdown




Training. Check this box if the event is a training session



Tap on Start •


The Incident Screen will appear.




Add Resources to an Open Incident

Choose the Assignment from the assignment dropdown or type an assignment into the


box



Scan Resources




Add a New Company to an Open Incident - Manual




Enter

Name*. Type in a company name which is required. The maximum number of characters is 50


ID*. Type in a company ID which is required. The maximum number of characters is 50.



The company will appear




Scan resources



Set Leader

Highlight the resource


Tap on the **Leader box**




The Leader will appear at the top of the list




Add a New Company to an Open Incident - Electronic




The company will appear in the incident and those resources that have been checked will appear





The Breaking Company box will appear




Remove Resources from a Company in an Open Incident - Electronic

Tap on the Off Scene assignment

Scan the resource



The Confirm Move screen will appear.




MOBILE ICS Functionality

The basis for the functionality of *interTRAX*TM Mobile is Incident Management System Accountability. Each Company and each Resource not attached to a Company <u>must be attached</u> to an Assignment. Each Evacuee or Patient <u>must be attached</u> to an Assignment. Assignments may be scanned or manually entered.


When you start a new incident the **Incident** screen will appear. The Mobile Assignment dropdown comes preloaded with the OFF SCENE and ON SCENE assignments

The Assignment that appears in the Assignment box is the active assignment and all companies, responders, evacuees or patients when entered will appear as part of that assignment.



Examples

- 1 Assignment
- 2 Scan
- 3 Information container
- 4 Resource attached to an assignment





















为不过不过其实是

Open, Close and End an Incident


Open

To open an incident tap on the incident to highlight and tap on the **Open** icon

An incident that is still in progress will have a YES status in the Open column.

An incident that has been ended will have a NO status in the Open column. You may open an ended incident.

NOTE: Only one incident may be open at a time.




Close or end on incident

To close or end an incident tap on OK in the bottom right hand corner of the screen. The End Incident box will appear.

Yes. The incident will be ended.

No. The incident will close and will not be ended.

Cancel. You will return to the open incident.



Using interTRAX™ MOBILE (EMS)

New Incident`

Tap on Menu

Choose New Medical Incident



The New Medical Incident screen will appear

Name. Type in the name or use the default Number. Type in the Number or use the default Incident. Choose from the dropdown. This is optional however it may be used to attach to an in progress interTRAX Exchange medical incident. Assignment. Choose from the dropdown. This is optional unless you will be attaching to interTRAX Exchange

Training. Check this box if the event is a training session



Tap on **Start**

The Mobile (EMS) Patient Log screen will appear

Detailed

Scan a triage tag or click on the No-Tag button enter patient information

Scan

Scan a triage tag. You may scan a driver license, assignment and/or equipment to add additional information or you may enter it manually.

Manual Entry

Tap on the **No-Tag** button to enter all of the information manually.



Enter the Patient information

The red items are required.



Tap on Save

The Patient Log screen will appear

The screen displays the Identity Code, Assignment and Triage Level

The patient remains highlighted in blue until the information has been synced.



Rapid Entry

Tap on the Rapid Entry button.




Enter the information.


The red items are required.

Scan any number of triage tags using the information that you entered.

When you tap on OK the information will be sent to the Patient Log.



To continue entering patients click on the **Begin** button.



Patient Details

Tap to highlight the Patient

Tap on the icon in the left hand corner

The Patient screen will appear.

You may change any of the information except the ID#

